

## Soli-Tone

### photo-energy system

A uniquely advanced technology for skin rejuvenation

#### Your treatment options:

LumiLift<sup>®</sup> treatment, ideal for lifting

- Face
- Eyes
- Chest

LumiFacial<sup>®</sup> treatments for face or body:

- Minor Acne
- Uneven pigmentation
- Redness
- Photo-aging

#### How long does a LumiLift<sup>®</sup> or LumiFacial<sup>®</sup> treatment take?

The regular Soli-Tone treatment takes only 30 minutes. However, if combined with a facial procedure and mask application the treatment will take one hour.

**Soli-Tone... visible results in one treatment!**



## Soli-Tone

a new multi-platform technology that reduces skin aging

#### LumiLift<sup>®</sup> treatment:

- Lessens fine lines and wrinkles
- Lifts and tones the facial contour



... all with no downtime or adverse effects!

LumiFacial<sup>®</sup> treatment uses a blend of wavelengths and different colored lights to correct skin conditions such as:

#### Acne



#### Brown age spot



#### Redness



#### Photo-damaged skin



Call now for your appointment

Available exclusively from your authorized trained esthetician

  
SILHOUET-TONE<sup>®</sup>

Silhouet-Tone 01 06 04 • code 420840



**SOLI-TONE**  
infinite energy  
**photo-energy system**

A uniquely advanced technology for skin rejuvenation



  
SILHOUET-TONE<sup>®</sup>



**The Soli-Tone is a unique skin rejuvenation device that uses light therapy and special currents to turn back the clock on aging, sun-damaged skin, minor acne, redness, and uneven pigmentation. This multi-platform system features two distinct facial treatments: the LumiLift® and the LumiFacial®. Simultaneous application of these technologies mutually enhances their individual effectiveness, promoting cellular rejuvenation from within. Results are immediate and visible.**



### **Are these treatments currently available?**

No. Up to now, treatments targeted individual problems. **Soli-Tone** is the first treatment device to combine two proven technologies: Light energy and microcurrents to treat conditions in the epidermis, dermis, and muscular regions without discomfort.

### **What is exactly the LumiLift® procedure?**

The **LumiLift®** procedure is a revolutionary way of rejuvenating the skin and lifting sagging muscles. It is performed with two (2) special electrodes that feature an advanced application of Micronized current and Light energy.

During the **LumiLift®** treatment, a pulsating light is directed at the skin surface. Light passes through the layers of the skin and penetrates into the dermis. This powerful light stimulates the fibroblast cell to produce more collagen. The increased production of collagen smoothes and softens the appearance of wrinkles, reduces the size of pores, and improves skin texture. At the same time, the electrodes emit Micronized currents to tone and lift sagging muscles.

At the end of the treatment, an antioxidant product is infused into your skin using two (2) special electrodes to smooth away facial creases.

Combining Light Therapy with Micronized currents erases years from the appearance of your skin without the use of drastic and painful surgical procedures. **LumiLift®** restores a natural, youthful appearance to the skin.

### **How does the LumiFacial® treatment work?**

The **LumiFacial®** targets common skin problems, such as minor acne, redness, and uneven pigmentation. During the treatment, a wrap-around facial panel is positioned in front of your face. The **Soli-Tone** selects the specific Light energy required for your skin type and condition.

At the end of the treatment, the esthetician will apply special products using a polarized current. This final step enhances the penetration of these products and promotes their effectiveness.

The treatment is completed after only 30 minutes. You can immediately return to your normal routine with no redness or peeling.

### **How does the treatment feel?**

The treatment is very comfortable; most clients fall asleep during the procedure. Unlike laser or IPL treatments that rely on thermal damage, **LumiLift®** and **LumiFacial®** treatments are non-invasive and totally natural. They use light from the visible spectrum, similar to sunlight but without the harmful UV rays.

### **How many treatments do I need?**

Age, lifestyle, and skin condition will influence the number of treatments needed. During the complimentary consultation, your skin type and condition will be determined. Your esthetician will be able to evaluate your skin, explain the needed procedures and tell you the estimated number of sessions needed.

While clients usually report a noticeable improvement after just one treatment, a series may be required. Generally, 6 to 10 treatments achieve optimal results.

### **What results can I expect?**

The improvement can be quite dramatic. Most clients see immediate improvement. With a series of treatment, years are easily removed from the appearance of the skin. You will see younger-looking skin, reduction in fine lines and wrinkles, increased skin elasticity and redefined and firmer facial contour. Pore size is improved along with textural irregularities, acne, and pigmentation blemishes.

This is achieved with no thermal damage to the skin and no down time. **Soli-Tone** procedures are unique in that the face can be helped as well as the neck, chest, and hands.

### **How long will the effects of the LumiLift® procedure last?**

Improvements continue for at least three (3) months after a single session. After completing a series, depending on the rate of the aging process, changes can last over one year. A good skin care regimen combined with a monthly maintenance treatment is sufficient to enhance and maintain the result of obtain with the **LumiLift®** treatment.

### **What are Micronized currents?**

Micronized currents are non-invasive. They stimulate and resemble the biological currents that your body produces to tighten muscles. Extensive clinical research has shown that these currents are able to lift and tone the muscles.

